




COMPACT4 Advanced line

Number of chambers	4
Max. pressure [mmHg]	220
Number of proprietary therapeutic pre-programs of the medical organization	6
Features for visually impaired/blind people	Voice assistant + special keyboard
Smart Gradient' function	mmHg, %, final pressure
Hold function [s]	1-50
Interval function [s]	1-50
Precise regulation of pressure and time	in steps of 1,5 or 10
Setting of individual pressure in each chamber	included
Number of independent channels	1
Adjustment of the speed of the treatment	1-5
Cuffs	overlapping system
User's own programs	included
List of "Favourites" programs	included
Display	5", TN LCD color touch screen
Weight [kg]	3,1
Dimensions [mm]	300x260x160
Power supply	100÷240 V, 50÷60 Hz
Power consumption	40-70 W
Warranty	1 year

Accessories




Included:




power supply

Optional:




4-chamber arm cuff

CPARM4T




4-chamber waist cuff

CPWAIST4T



4-chamber left leg cuff

CPLEG4TL




4-chamber right leg cuff

CPLEG4TR




4-chamber short pants

CPSHORTPANTS4T




Extension zipper

CPLEGEXT4




Double splitter

CPDOUBLECONNECTOR4




Patient's safety button

SAFETY BUTTON



Double splitter (long - 150 cm)
for 4-chamber devices

CPDOUBLECONNECTOR4L



2 x Extension zippers
for 4-chamber short pants

CPSHORTPANTSEXT4

Therapeutic pre-programs

of the medical organization



lymphoedema |

-
-
-
-

Gradient	1 mmHg
Time	50 min.
Hold	3 s
Interval	3 s

venous oedema |

-
-
-
-

Gradient	3 mmHg
Time	45 min.
Hold	2 s
Interval	3 s

lipoedema |

-
-
-
-

Gradient	2 mmHg
Time	60 min.
Hold	3 s
Interval	3 s

regeneration

-
-
-
-

Gradient	1 mmHg
Time	60 min.
Hold	3 s
Interval	3 s

prevention, relaxation

-
-
-
-

Gradient	0 mmHg
Time	60 min.
Hold	2 s
Interval	2 s

all-chamber massage

-
-
-
-

Gradient	0 mmHg
Time	30 min.
Hold	2 s
Interval	3 s